

# Proteomics Using Mass Spectroscopy

## Can the FEL Help?

Dariya Malyarenko

Dennis Manos


Michael Trosset

Gene Tracy


Haijian Chen

Pete Harris

Christine Hopkins





# Surface Enhanced Laser Desorption/Ionization (SELDI) Time-Of-Flight (TOF) Mass Spectroscopy



# W&M Mass Spectroscopy Proteomics Goals

- Enhance Signal/Noise
- Connect markers to biological processes



# Data processing




- Resampling improves peak identification

# Can the FEL help?


- Analog electronics and detectors
- High rep rate FEL could lead to ion counting



# Haglund et al. use the Vanderbilt FEL for MALDI



- TOP: IR-MALDI spectrum of cytochrome-C in succinic acid (single shot). Note the appearance in the spectrum of multiply-charged ions with high efficiency.




- BOTTOM: Ten-shot MALDI spectrum of insulin in succinic acid using a special fast-drying matrix preparation to make micron-size crystallites. Note the absence of matrix background. Mass resolution  $m/Dm \sim 600$ .

# Ionization is complicated!


- Large pulse energies
- Collisional ionization
- Big spots needed!

B.J. Garrison, A. Delcorte, L.V. Zhigilei, T.E. Itina, K.D. Krantzman, Y.G. Yingling, C.M. McQuaw, E. J. Smiley, and N. Winograd, *Appl. Surf. Sci.* 203-204, (2003).



# Can the FEL help?


- Source photoionization to enhance the signal
- High rep rate and high power could make single ion counting possible (everything ionizes at a  $10^{14}$  W/cm<sup>2</sup>)
- What is the role of wavelength?






# Limits to SNR


- Detector gain variation




- FEL could fragment/ionize at detector




e production depends  
on velocity, not  
energy



# Will the FEL enable higher resolution?



(c)




(e)

- Filtering resolves hidden peaks
- Can we physically separate them?

# Shifted peaks can mean different biology


- Proteins come with modifications:
  - Phosphorylated
  - Acetylated
  - Methylated
  - Sulfated




How can you tell if a shoulder is a peak or a problem?

# Summary of Biomarker Discovery and Identification

## SELDI-TOF




## Classification and Regression Tree Analysis




[Cancer Research 62, 3609-3614, July 1, 2002] © 2002  
[American Association for Cancer Research](http://www.aacr.org/)


## Identification



## Purification




# MS/MS for protein identification



# MS/MS for Using FEL Fragmentation/Ionization

With FEL  
fragmentation/ionization you  
can get more than one ion per  
parent!

- Timing determines the parent mass
- Vertical position determines the fragment mass



# The upgraded FEL is the perfect post-ionizer

- High rep rate makes it quasi continuous

$$\Delta x = (30000 \text{ m/s})(12 \text{ ns}) = 360 \mu\text{m}$$

As long as the Rayleigh range is at least 360 microns, every ion will see a pulse.

# The upgraded FEL is the perfect post-ionizer

- High peak power enhances multi-photon ionization or fragmentation

$$I_{\text{max}} = \frac{2 \text{ kW}}{10^{-4} \text{ cm}^2} \left( \frac{12 \text{ ns}}{0.2 \text{ ps}} \right) \times 100 = 10^{14} \text{ W/cm}^2$$



# The upgraded FEL is the perfect post-ionizer

- Tunability or intensity variation may enhance biological selectivity
- Specifics of fragmentation are unimportant, as long as there is a difference!

# Conclusion

- High throughput MS/MS will be possible using the upgraded FEL.
- MS/MS will simultaneously improve resolution and protein identification.
- CW laser post-ionization requires understanding of the ionization and fragmentation processes for biomolecules
- New fast, imaging detectors are necessary